

AURORA **HD**

My dream to see
the **vibrant world**
again

Offering superior picture quality and a large viewing area, Aurora HD is ideal for reading larger books, writing or doing hobbies.



Plug and play



Movable monitor



High definition image




ZOOMAX[®]
www.zoomax.co

Technical Specifications

Display	24" LCD, 1920 x 1080 pixels
Camera	High definition and auto focus, active picture elements of the camera: 1280 x 720pixels
Magnification	2X ~ 70X
Color modes	17 alternative color modes 5 default color modes
Basic functions	Zoom in/out, Color modes adjustment Brightness adjustment
Advanced functions	Freeze frame, Focus lock, Find feature Reading lines and windows, Memory setting
Workspace under the camera	18cm
Ergonomic adjustment	Customized position of light and adjustable illumination
Weight	15.8kg
Certificates	CE; FCC; RoHS



Contact us

 Tel: +86-571-87006308
 Email: sales@zoomax.co
 Web: www.zoomax.co