

Panda

See simpler live easier

With a compact design and competitive price, Panda is perfect for a home or apartment with limited space and any budget.



 3 buttons only

 Space-saving size

 High quality with affordability


ZOOMAX[®]
www.zoomax.co

Technical Specifications

Display	19" LCD
Camera	Auto focus
Magnification	2x ~ 70x
Color modes	17 modes
Favorite color modes	Full color + 4 favorite color modes
Ergonomic adjustment	Screen height adjustable
Functions	Zoom in/out, switch color mode Adjust brightness, save last setting
Working height	25 cm
XY table dimensions	36.8cm×36.6cm×4.6cm X dimension travel: 20cm Y dimension travel: 23cm
Weight	17.6kg
Certificates	CE; FCC; RoHS



Contact us

 Tel: +86-571-87006308
Email: sales@zoomax.co
Web: www.zoomax.co