




NEW

Snow 10 PRO

Feel More Beyond Big

Zoomax Snow 10 Pro is the perfect blend of portability, affordability, and advanced functionality typically found in devices costing much more. Standard features include OCR and speech, and a 3-way camera that adjusts for distance viewing, reading, or as a mirror. With an MP3 music player, the ability to store 1,000 images, and magnification to 19x, Snow 10 Pro delivers!



-  20 languages text-to-speech
-  Reader mode for easier reading
-  Self-view/distance view/close view

Key Specifications

Screen	10" touch screen; screen resolution: 1280 x 720p
HD Camera	13M pixels; auto focus/tap to focus
Color modes	Full color & 10 high contrast color modes
Viewing modes	Self/Close/Distance
OCR & Speech	Optional-20 languages
Magnification	2x ~ 19x (non-text-to-speech) 1.3x ~ 8x (text-to-speech) 1.3x (reader mode)
Storage	1.5GB (up to 1000 images)
Port	Type-C/Micro-HDMI/Headphone
Fast charging time	2 hours
Battery	Approximately 3.5 hours of continuous use
Weight	2 lbs/937g/33oz
Dimensions	280 x 190 x 20.9mm (11 x 7.5 x 0.82in)
Certificates	CE FC RoHS

