

# A DAY WITH SNOW 7 HD PLUS

Just another peaceful day.

## MORNING

Fully Charged for another great day!



## TAKE A BUS

Never saw the bus line so clear before.



## OFFICE

Now I can do the paperwork in high contrast color!



## CHAT TIME

What's happening around today?



## SHOPPING

What a big discount! I must try this!



## DINNER

The Speech function can read my new recipe out.



## READING TIME

Good, I can use EaseArm for speech and writing!







## GOOD NIGHT


Before bed, I can transfer images from PC to Snow 7 HD Plus.




## SNOW 7 HD PLUS SPECIFICATION

Magnification	2.5X ~ 19X	
Display Modes	Full color and 10 high contrast color modes	
Reading Line	Fixed horizontal/vertical reading line	
Screen	7" touch screen	
Screen Resolution	1280x720p	
Camera	13M pixels	
Focus	Auto focus/tap to focus	
Viewing Distance	5 m	
Connection	Connection to TV / Monitor through HDMI cable	
Image Transfer (through USB cable)	Transfer from PC/Mac to device Transfer from device to PC/Mac	
Storage	For up to 1,000 photos	
Battery Charging Time	Approximately 3.5 hours 3.5 hours (+/-15 minutes)	
Standby Mode	Short press the Power button or stop using after 3 minutes	

Weight	532 g / 18.8 oz.
Dimensions	195 mm x 139.5 mm x 25 mm / 7.7" x 5.5" x 0.98"

<b>Text-to-Speech</b>		
Magnification	1.7X ~ 8X 7 magnification levels in the text page	
Language	20 (more are coming)	
Voice Gender	Female/Male (6 languages available only in female)	
Full Page Reading	Available	
Text Color	4 high contrast color modes	
EaseArm (optional)	Swivel arm facilitating writing and text-to-speech	

<b>Power Adapter</b>		
Input	100 - 240V	
Output	5V/3.5A	

<b>Operation and Storage Condition</b>		
Humidity	< 70%	
Temperature	10-40 °C	