


# Snow 4.3"


The World is Clear Again

Zoomax Snow 4.3 inch handheld video magnifier is given advanced image technology and camera. It offers sharp image even at the lowest magnification and under fast movement. With large buttons and the unique one hand control design, it's very easy to use even without any learning at the very beginning. You can read and write on the go or at home.



 Easy to use with intuitive buttons

 Light weight

 Big and intuitive buttons

## Technical Specifications

Display	4.3" TFT screen
Camera	Auto focus
Magnification	2x to 16x
Color modes	10 modes
Image freeze	Yes
Writing stand	Yes
Connect to TV	PAL / NTSC
Battery	3.5 hours Li-ion battery, replaceable & rechargeable
Power saving mode	Enter into power saving mode automatically when it has not been used for over 3 minutes
Weight	204g (7.20oz) with battery / 165g (5.82oz) without battery
Dimensions	160mm x 80mm x 24.8mm (6.30" x 3.15" x 0.98")
Certificates	CE FC RoHS

