



Cost-effective Device?

Widescreen?

See Near and Far?

**SNOW 7 HD**, the first 7" video magnifier you feel really easy to use!



◀ Near Viewing



◀ Distance Viewing

Folding Frequency of Reading Stand

**8000 Times**

Tailored Safety Battery with

**Temperature Protection**

# Lightest 7" Video Magnifier



**2 Mins** Learning to Use

**480g**



**Crisp Image Instantly**



**Ergonomic Reading Stand**

Reduce neck strain



**Yes!**


SNOW 7 HD



**No!**

Other 7"

**5 Large Buttons Only**



**Battery Life**

Approximately 5 hours of continuous use



**Special Protective Bag**



## Technical Specifications

### Magnification

- 2.2X ~ 19X

### High Contrast Color Modes

- 12 modes (3 favorite color modes)

### Screen

- 7-inch TFT

### Camera

- Auto focus: 5M pixels

### Viewing Distance

- 6.5cm ~ 500cm

### Functions

- Freeze image; Switch color mode; Zoom in/out; Setting memory

### Save Power

- The Snow 7 HD will turn off after 3 minutes of non-use

### Connect to TV/Monitor

### Power Adapter

- Input: 110-240V
- Output: 5V/2.5A; Power: 2.5A

### Li-ion Battery Life

- Approximately 5 hours

### Charging Time

- 4~5 hours

### Dimension

- 19.5cm×13.95cm×2.5cm

**ZOOMAX**<sup>®</sup>  
Zoom Your Vision, Max Your Life



Tel: +86-571-87006308  
Email: sales@zoomax.co  
Web: www.zoomax.co

**SNOW7 HD**  
7" HD Handheld Video Magnifier



*Easy is  
what you want!*

ZOOMAX, a Rejoin Solutions' subsidiary, is a developer and manufacturer of technology dedicated to creating independence for people with vision impairments, offering a full line of handheld, portable and desktop video magnifiers.

**ZOOMAX**<sup>®</sup>  
www.zoomax.co