

Tight Budget?

Magnifier Too Heavy?

Want Crisp Images?

Seek Portability?

Want to See Near and Far?

Desire Multi-connectivity?

Mars HD, offers reading, distance and self-viewing modes in one compact device. And all at an unbeatable price!



◀ School



◀ Office



◀ Home

Lightest Portable Video Magnifier



Set-up in 5 Seconds

Less than 0.9 kg



Small Size
Perfect for any environment



Crisp Image Instantly



Battery Life
4 hours of continuous use



Convenient Carrying Case



Optional

Multi-connectivity

Connect to Mac, PC or monitor through USB and HDMI port



Technical Specifications

Camera

- 12X optical zoom, Autofocus
- HD with resolution 1280 x 720 pixels
- Rotate 330° vertically or 180° horizontally for multi-angles reading

Frame Rate

- 60 full digital frame /second @ 720p

Working Height Under Camera

- 19.5cm

Magnification Ranges

- 2X ~ 80X (based on 14" screen)

Color Modes

- Full color + 16 high contrast color modes

3 Viewing Modes

- Distance viewing; Self viewing; Reading viewing (left/right)

USB/HDMI Connection

- Connect to Mac, PC or monitor

Functions

- Zoom in/out, Color mode adjustment, Brightness adjustment, Reading lines/windows, Freeze frame, Split screen, Find, Image saving, Led light on/off, Setting memory

Li-ion Battery

- Rechargeable, 4 hours continuous use

Net Weight

- Less than 0.9kg

ZOOMAX[®]
Zoom Your Vision, Max Your Life



Tel: +86-571-87006308

Email: sales@zoomax.co

Web: www.zoomax.co

MARS HD
Portable Video Magnifier



ZOOMAX, a Rejoin Solutions' subsidiary, is a developer and manufacturer of technology dedicated to creating independence for people with vision impairments, offering a full line of handheld, portable and desktop video magnifiers.

ZOOMAX[®]
www.zoomax.co